

OmniMax™

MMF System



CMF



ZIMMER BIOMET
Your progress. Our promise.™

We reinvented the bar. Then we raised it.

Bone-borne arch bars offer the strength and stability of Erich arch bars combined with the speed and simplicity of IMF screws. They also provide surgeons, hospital staff and patients with several benefits:

- Reduced need for interdental wiring
- Minimized risk for penetrating wire-stick injury
- Increased OR efficiency may occur due to shortened application/removal time
- Expanded application and removal options – procedure may be performed in a hospital or office setting

It's not just what we make,
**It's what
we make possible**



1 Extended Screw Insertion Slots


Our elongated, semi-locking slotted arch bar allows for variability in screw placement for tooth root avoidance.

2 Type II Anodization

OmniMax Arch Bars are Type II anodized to increase fatigue strength and offer a smooth surface that minimizes sharp edges which may reduce patient irritation.

3 Adjustable, Semi-Locking Technology

The unique arch bar and screw engage together to elevate the arch bar off of the soft tissue and to help prevent gingival compression.



4 An All-Inclusive, Adaptable MMF System


- The OmniMax MMF System includes bone-borne arch bars, Erich Arch Bars and uniquely-designed IMF screws.
- The system can be used as a stand-alone tray and is also compatible with the TraumaOne™ system for MMF procedures.
- The outer container can be easily customized to house additional surgical instrumentation or implants - adapting to your surgical preferences.

Step 1. Determine Arch Bar Length



Determine length of arch bar required and cut to appropriate length, if required.

TIP: A deburr area can be found on the side of the Adjustment Tool, item #01-0292, for smoothing of rough edges after cutting.



Deburr area on adjustment tool.


Step 2. Contour The Arch Bar



Fig. 1



Fig. 2



Manually contour the arch bar to approximate the shape of the maxilla or the mandible.


TIP: The arch bar may require additional contouring in the mandible to accommodate patient anatomy. See image for guidance. Perform slight in-plane bend to most posterior segment on each side (**Fig. 1**). Perform slight twist to most posterior segment on each side (**Fig. 2**).

Step 3. Insert Screws




Insert the self-drilling screw into the bone through a single slot of the arch bar, in between tooth roots, until the bar is seated in the retention groove.



Tip: For screw insertion at angles in excess of 10 degrees, the Adjustment Tool helps to enable proper engagement into the bar.



Step 4. Adjust Gingival Height



As screw engages with the arch bar, the bar is pulled towards the head of the screw until seated in the retention groove (Fig. 1, Fig. 2). Upon seating in the retention groove, the final soft tissue stand off can be adjusted by inserting the screw closer to or further away from the gingiva (Fig. 3, Fig. 4).



Step 5. Repeat Screw Insertion



Repeat screw insertion along the arch bars, ensuring that a minimum of three screws are used per bar and no more than two consecutive slots remain empty. Four screws per bar may be required depending on the length of the arch bar used.

Step 6. Wire Bars Together



A minimum of two, 24-gauge ligature wires should be wrapped around the arch bar hooks between the maxilla and mandible to achieve MMF. Different gauges (up to 18 gauge wire) of wires or elastics may be used to achieve MMF. **Arch Bar Removal:** The OmniMax MMF System can be removed in the office or clinic with local anesthesia or local anesthesia with sedation at 4 to 6 weeks postoperative.

Implants and Instruments

2.0mm OmniMax and IMF Screws

Self-drilling OmniMax MMF Screws






Part No.	Description
91-5707	2.0 x 7.0mm
91-5709	2.0 x 9.0mm
91-5711	2.0 x 11.0mm

Self-drilling IMF Screws



Part No.	Description
91-5607	2.0 x 7.0mm
91-5609	2.0 x 9.0mm
91-5611	2.0 x 11.0mm

Bars and Wires



Blades and Drills



Instruments



OmniMax MMF Adjustment Tool
with Deburr
01-0292



Wire Cutter
51-0928



Wire Twister
51-6705



U of M Tongue Depressor
and Cheek Retractor
09-0103*



Wieder Retractor
09-0105*



Black Ratcheting Handle
46-0008*

Containers and Tray Components



OmniMax MMF Tray
46-2462



OmniMax MMF Outer Container
46-2463



OmniMax MMF Auxiliary Tray
46-2464



The OmniMax MMF System is compatible with the TraumaOne system for MMF procedures.



ZIMMER BIOMET

Your progress. Our promise.™

For more information on OmniMax™ MMF and other craniomaxillofacial solutions, please contact us at:

BIOMET MICROFIXATION GLOBAL HEADQUARTERS

1520 Tradeport Drive • Jacksonville, FL 32218-2480

Tel 904.741.4400 • Toll-Free 800.874.7711 • Fax 904.741.4500 • Order Fax 904.741.3059

www.zimmerbiomet.com

EUROPE

Toermalijnring 600 • 3316 LC Dordrecht • The Netherlands

Tel +31 78 629 29 10 • Fax +31 78 629 29 12

All content herein is protected by copyright, trademarks and other intellectual property rights owned by or licensed to Biomet, Inc. or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet. This material is intended for healthcare professionals and the Zimmer Biomet sales force. Distribution to any other recipient is prohibited. Zimmer Biomet does not practice medicine. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this information does not replace the need for such examination and/or advice in whole or in part. Each physician should exercise his or her own independent judgment in the diagnosis and treatment of an individual patient. For complete product information, including indications, contraindications, warnings, precautions and potential adverse effects, see the package insert and www.zimmerbiomet.com.