



RibFix Blu

THORACIC FIXATION SYSTEM


Closing the Gap in Rigid Fixation

RibFix Blu is a comprehensive thoracic fixation system that allows for the stabilization and rigid fixation of fractures in the chest wall including sternal reconstructive surgical procedures, trauma or planned osteotomies. The system may be used in normal and poor bone to promote bone union.

Features and Benefits: RibFix Blu Thoracic Fixation System features implants and instrumentation designed to improve OR efficiency.

Innovative plate-to-bone approximation tools

- Reduce plates flush to bone
- Allows for precise placement of plates along the rib



Self-drilling screws

- Eliminates need for pre-drilling
- Improves OR efficiency





MIS trans-scapular instrumentation

- Standard options for percutaneous, small incision, sub or trans-scapular access to hard-to-reach fractures

Unique plate contouring tools

- Make adjustments to the plate in all planes, including torsional bends
- Eliminates need to remove implant from surgical field



RibFix Blu Thoracic Fixation System

