





DORO LUCENT[®] Disposable Skull Pins





DORO LUCENT® Disposable Skull Pins



Titanium Skull Pins

DORO LUCENT®

Disposable Skull Pins

Item no. 1106.003 (3 pins)

Item no. 1106.103 (9 pins)

- For intra-operative imaging applications using CT, MRI, X-Ray and Fluoroscopy.
- Designed to reduce imaging artifacts to an absolute minimum.
- The composition of titanium tip and high-performance reinforced polymer guarantees a precise, secure and rigid fixation.
- Sterile packaging.

Headquarters Germany

pro med instruments GmbH | Bötzingen Straße 38 | 79111 Freiburg, Germany
Phone +49 761 384 222 10 | E-mail pmi@pmisurgical.com

Version 02 © 2017 pro med instruments GmbH, Freiburg, Germany. All rights reserved.

The information has been compiled with utmost care. Errors can, however, not be excluded. Errors do not entitle to any legal claims.
We reserve the right to modify the products at any time.

US Subsidiary

pro med instruments, Inc. | 4529 S.E. 16th Place, Suite 101 | Cape Coral, FL 33904
Toll Free 877 225 4086 | E-mail pmi.us@pmisurgical.com