

iQ[®] Intelligent System

Neurosurgery



CMF



**You've concentrated on the procedure.
Let it concentrate on the closing.**

Unprecedented Speed


The iQ inserts screws in less than a fraction of a second. It's the fastest insertion method available to neurosurgeons today.

Automated Ease

A slight downward force and the push of a button is all it takes to insert screws with precision and optimal torque every time.

Intelligent Assistance

The iQ's internal software carefully monitors torque and controls the motor to ensure consistent screw insertion and prevent screw stripping.*



Experience faster screw insertion with the iQ Intelligent System.*

- When compared to manual screw insertion, the iQ provided time savings triple that of battery-powered drivers.
- Screw insertion using the iQ is 40% faster than traditional battery-powered screwdrivers.
- Screw insertion using the iQ may reduce time and cost by contributing to a faster closure.

*Evaluation of the iQ Intelligent System for Rapid Screw Insertion, Hatcher, B

To learn more about the iQ Intelligent System, contact your Zimmer Biomet representative or call 1.800.874.7711.